

# PHARMA & LAB REFRIGERATOR



## USP

- AIR CIRCULATING FAN
- AUTO DEFROST SYSTEM
- ERGONOMIC HANDLE WITH LOCK & KEY
- FLAT & STURDY SHELVES
- DIGITAL MICRO CONTROLLER PROCESSOR WITH LED DISPLAY
- 4 WHEELS WITH BRAKE

## Laboratory Refrigerator Double Door (2°C To 8°C)

- Air Circulating Fan
- Auto Defrost System
- Ergonomic Handle with Lock & Key
- Flat 7 Sturdy Shelves
- Digital Micro Controller Processor with LED Display
- 4 Wheel with Brake
- High/low temperature audio-visual alarm system for safety

Item Code	Capacity (Ltr.)	External (Dimension)	Internal (Dimension)	No. of Shelves
JW-1000	1000	1025 X 745 X 2030	941 X 648 X 1581	14