

PHARMA & LAB REFRIGERATOR




Laboratory Deep Freezer (-40°C)

- Build for tropical climates, tested at 41°C / 75% RH ambient
- Controlled reliable functionality with fitted in controller
- Wide display area
- Low maintenance condenser
- Easy serviceability with quick electrical connectors and easy access components
- Fast defrost in and around 2 minutes

Item Code	Capacity (Ltr.)	External (Dimension)	Internal (Dimension)	No. of Shelves
-----------	-----------------	----------------------	----------------------	----------------

JW-IDF-300	300	576 X 650 X 1454	492 X 551 X 1050	4
JW-IDF-400	400	672 X 690 X 2046	552 X 535 X 1400	5
JW-IDF-500	500	672 X 690 X 2200	552 X 535 X 1660	5

* Vertical Deep Freezer having cabin insulation 60mm with Door Lock.



- Tropical ideal for 43°C
- IP24
- All Components Safety Certified
- European design and Engineering
- Corrosion resistant body
- High energy efficiency compressor
- Anti rodent protection system
- Maintenance Free Condenser
- Soft look and robust design