

PHARMA & LAB REFRIGERATOR

Laboratory Refrigerator Single Door (2°C To 8°C)

- Air Circulating Fan
- Auto Defrost System
- Ergonomic Handle with Lock & Key
- Flat 7 Sturdy Shelves
- Digital Micro Controller Processor with LED Display
- 4 Wheel with Brake
- High/low temperature audio-visual alarm system for safety

Item Code	Capacity (Ltr.)	External (Dimension)	Internal (Dimension)	No. of Shelves
RE-200	200	450 X 539 x 1464	380 X 425 X 1100	5
RE-250	250	450 X 539 X1724	380 X 425 X 1350	5
RE-300	300	576 X 650 X 1454	492 X 551 X 1050	7
RE-400	400	576 X 650 X 1724	492 X 551 X 1320	7
RE-450	450	576 X 650 X 2055	492 X 551 X 1664	7
RE-500	500	620 X 676 X 1897	540 X 577 X 1660	8



Laboratory Deep Freezer (-20°C)

- Build for tropical climates, tested at 41°C / 75% RH ambient
- Controlled reliable functionality with fitted in controller
- Wide display area
- Low maintenance condenser
- Easy serviceability with quick electrical connectors and easy access components
- Fast defrost in and around 2 minutes

Item Code	Capacity (Ltr.)	External (Dimension)	Internal (Dimension)	No. of Shelves
RE-DF-300	300	576 X 650 X 1454	492 X 551 X 1050	4
RE-DF-400	400	672 X 690 X 2046	552 X 535 X 1400	5
RE-DF-500	500	672 X 690 X 2200	552 X 535 X 1660	5

* Vertical Deep Freezer having cabin insulation 60mm with Door Lock.